

## IHE INTEGRATION STATEMENT

Vendor	Product Name	Version	Date
Candelis, Inc.	ImageGrid	1.7	April 2010

This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:

Integration Profiles Implemented	Actors Implemented	Options Implemented
Access to Radiology Information	Image Manager/ Image Archive	
	Report Repository	
Consistent Presentation of Images	Image Manager/ Image Archive	
Evidence Document	Image Manager/ Image Archive	
Import Reconciliation Workflow	Image Manager/ Image Archive	
	PPS Manager	
	DSS/OF	
Key Image Note	Image Manager/ Image Archive	
Mammography Image	Image Manager/ Image Archive	
Scheduled Workflow	Image Manager/ Image Archive	
	PPS Manager	
	DSS/OF	Assisted Acquisition Protocol Setting Option Modality Billing Material Option