

## IHE INTEGRATION STATEMENT

Vendor	Product Name	Version	Date
Candelis, Inc.	ImageGrid Mammography Viewer	1.7	April 2010

This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:

Integration Profiles Implemented	Actors Implemented	Options Implemented
Schedule Workflow	Image Display	No options defined
Mammography Image	Image Display	None
Access to Radiology Information	Image Display	None

**Internet address for vendor's IHE information:** <http://www.candelis.com/ihe>

### Links to Standards Conformance Statements for the Implementation

<b>HL7</b>	<a href="http://www.candelis.com">http://www.candelis.com</a>	
<b>DICOM</b>	<a href="http://www.candelis.com">http://www.candelis.com</a>	

### Links to general information on IHE

In North America : <a href="http://www.rsna.org/IHE">http://www.rsna.org/IHE</a>	In Europe: <a href="http://www.ihe-europe.net">http://www.ihe-europe.net</a>	In Japan: <a href="http://www.jira-net.or.jp/ihe-j">http://www.jira-net.or.jp/ihe-j</a>
---	---	--