

<b>IHE Integration Statement</b>			
<b>Vendor</b>	<b>Product Name</b>	<b>Version</b>	<b>Date</b>
Candelis Inc.	ImageGrid	2.0	1/26/13

This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below :

<b>Integration Profile</b>	<b>Actor</b>	<b>Integration Profile Option</b>
Access to Radiology Information	Image Manager/Archive	None
Audit Trail and Node Authentication	Secure Node	None
Consistent Presentation of Images	Image Manager/Archive	None
Consistent Time	Time Client	None
Evidence Documents	Image Manager/Archive	None
Key Image Note	Image Manager/Archive	None
Mammography Image	Image Manager/Archive	None
Portable Data for Imaging	Portable Media Creator	DVD Media
Portable Data for Imaging	Portable Media Creator	None
Patient Information Reconciliation	Image Manager/Archive	None
Patient Information Reconciliation	Department System Scheduler/Order Filler	None
Patient Information Reconciliation	Performed Procedure Step Manager	None
Patient Identifier Cross-referencing for MPI	Patient Identifier Cross-reference Consumer	None
Patient Identifier Cross-referencing for MPI	Patient Identity Source	None
Scheduled Workflow	Image Manager/Archive	PPS Exception Management
Scheduled Workflow	Image Manager/Archive	Performed Work Status Update Receive

Internet address for vendors IHE information : <a href="http://www.candelis.com">http://www.candelis.com</a>	
Links to Standards Conformance Statements for the Implementation	
<b>HL7</b>	<a href="http://www.candelis.com/admin/uploads/ImageGrid_HL7ConformanceStatement.pdf">http://www.candelis.com/admin/uploads/ImageGrid_HL7ConformanceStatement.pdf</a>
<b>DICOM</b>	<a href="http://www.candelis.com/admin/uploads/Candelis%20DICOM%20Conformance%20Statement%202013.pdf">http://www.candelis.com/admin/uploads/Candelis%20DICOM%20Conformance%20Statement%202013.pdf</a>
Link to general information on IHE : <a href="http://www.ihe.net">www.ihe.net</a>	



IHE Integration Statement			
Vendor	Product Name	Version	Date
Candelis Inc.	ImageGrid	2.0	1/26/13

This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below :

Integration Profile	Actor	Integration Profile Option
Scheduled Workflow	Image Manager/Archive	Availability of PPS-Referenced Instances
Scheduled Workflow	Image Manager/Archive	None
Scheduled Workflow	Department System Scheduler/Order Filler	PPS Exception Management
Scheduled Workflow	Department System Scheduler/Order Filler	Image Availability Query
Scheduled Workflow	Department System Scheduler/Order Filler	Performed Work Status Update Receive
Scheduled Workflow	Department System Scheduler/Order Filler	Departmental Appointment Notification
Scheduled Workflow	Department System Scheduler/Order Filler	Availability of PPS-Referenced Instances
Scheduled Workflow	Department System Scheduler/Order Filler	None
Scheduled Workflow	Performed Procedure Step Manager	None

Internet address for vendors IHE information : <a href="http://www.candelis.com">http://www.candelis.com</a>	
Links to Standards Conformance Statements for the Implementation	
<b>HL7</b>	<a href="http://www.candelis.com/admin/uploads/ImageGrid_HL7ConformanceStatement.pdf">http://www.candelis.com/admin/uploads/ImageGrid_HL7ConformanceStatement.pdf</a>
<b>DICOM</b>	<a href="http://www.candelis.com/admin/uploads/Candelis%20DICOM%20Conformance%20Statement%202013.pdf">http://www.candelis.com/admin/uploads/Candelis%20DICOM%20Conformance%20Statement%202013.pdf</a>
Link to general information on IHE : <a href="http://www.ihe.net">www.ihe.net</a>	