

IHE Integration Statement

Vendor	Product Name	Version	Date
Candelis Inc.	ImageGrid	2.0	January 26,

This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:

Profiles Implemented	Actors Implemented	Options Implemented
Access to Radiology Information	Image Manager/Archive	
Consistent Presentation of Images	Image Manager/Archive	
Consistent Time	Time Client	
Evidence Documents	Image Manager/Archive	
Key Image Note	Image Manager/Archive	
Patient Information Reconciliation	Department System Scheduler/Order Filler	
	Order Placer	
	Performed Procedure Step Manager	
Portable Data for Imaging	Portable Media Creator	DVD Media
Scheduled Workflow	Department System Scheduler/Order Filler	
	Image Manager/Archive	
	Order Placer	
	Performed Procedure Step Manager	

Internet address for vendor's IHE information:

Links to Standards Conformance Statements for the Implementation

HL7	http://www.candelis.com/admin/uploads/pdf_files/imagegrid-hl7-1.4.72.pdf
DICOM	http://www.candelis.com/admin/uploads/pdf_files/imagegrid-dicom-1.4.72

Links to general information on IHE

In North America: www.ihe.net	In Europe: www.ihe-europe.net	In Japan: www.ihe-j.org
--	---	--