



| IHE Integration Statement | | | |
|--|--|--|--------------|
| Vendor | System Name | System Version | Date |
| CANDELIS | ImageGrid | 1.9 | 2012, Jan 26 |
| This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below: | | | |
| Integration Profile implemented | Actor implemented | Integration Profile Option implemented | |
| Access to Radiology Information | Image Manager/Archive | None | |
| Access to Radiology Information | Report Repository | None | |
| Consistent Presentation of Images | Image Manager/Archive | None | |
| Consistent Time | Time Client | None | |
| Evidence Documents | Image Manager/Archive | None | |
| Key Image Note | Image Manager/Archive | None | |
| Mammography Image | Image Manager/Archive | None | |
| Portable Data for Imaging | Portable Media Creator | DVD Media | |
| Portable Data for Imaging | Portable Media Creator | None | |
| Patient Information Reconciliation | Image Manager/Archive | None | |
| Patient Information Reconciliation | Department System Scheduler/Order Filler | None | |
| Patient Information Reconciliation | Order Placer | None | |
| Patient Information Reconciliation | Performed Procedure Step Manager | None | |
| Simple Image and Numeric Report | Enterprise Report Repository | None | |
| Simple Image and Numeric Report | Report Repository | None | |
| Scheduled Workflow | Image Manager/Archive | None | |
| Scheduled Workflow | Image Manager/Archive | Availability of PPS-Referenced Instances | |
| Scheduled Workflow | Image Manager/Archive | PPS Exception Management | |
| Scheduled Workflow | Image Manager/Archive | Performed Work Status Update Recieve | |
| Scheduled Workflow | Department System Scheduler/Order Filler | PPS Exception Management | |
| Scheduled Workflow | Department System Scheduler/Order Filler | None | |
| Scheduled Workflow | Department System Scheduler/Order Filler | Image Availability Query | |
| Scheduled Workflow | Department System Scheduler/Order Filler | Performed Work Status Update Recieve | |
| Scheduled Workflow | Department System Scheduler/Order Filler | Departmental Appointment Notification | |

| | | |
|---|---|--|
| Scheduled Workflow | Department System Scheduler/Order Filler | Availability of PPS-Referenced Instances |
| Scheduled Workflow | Order Placer | None |
| Scheduled Workflow | Order Placer | Departmental Appointment Notification |
| Scheduled Workflow | Performed Procedure Step Manager | None |
| | | |
| Internet address for vendors IHE information: http://www.candelis.com | | |
| | | |
| Links to Standards Conformance Statements for the Implementation | | |
| HL7 | http://www.candelis.com/admin/uploads/ImageGrid_HL7ConformanceStatement.pdf | |
| DICOM | http://www.candelis.com/admin/uploads/ConformanceStatement2010.pdf | |
| Links to general information on IHE | | |
| In North America: www.ihe.net | In Europe: www.ihe-europe.org | In Japan: www.jira-net.or.jp/ihe-j |