

StreamVue™ Enterprise Radiology Viewer

INFO SHEET

StreamVue is an enterprise radiology viewer engineered for instant access, uncompromising performance, and seamless integration.


StreamVue's advanced streaming technology delivers real-time, high-fidelity images to any device, anywhere, with no compromise on quality or security. Images load instantly, enabling radiologists to move through studies up to 30% faster and dedicate more time to diagnostic accuracy and patient care.

StreamVue is designed to work with the infrastructure you already have. Whether your PACS is deployed on-premise or in the cloud, StreamVue is compatible with any third-party system via standard DICOM, delivering a modern viewing experience without requiring a costly and disruptive PACS replacement.

By eliminating the need for specialized workstations, StreamVue extends diagnostic-quality access to radiologists wherever they work — in the reading room, at a satellite facility, or remotely — through either a full-featured multi-monitor application or a zero-footprint browser-based viewer.

Ideal for both small practices and large health systems, StreamVue is a complete diagnostic solution, equipped with the comprehensive suite of viewing, measurement, and workflow tools healthcare professionals expect from an enterprise-grade viewer.

StreamVue also supports the full spectrum of medical imaging modalities across vendors, ensuring clinicians have uninterrupted access to all imaging data while maintaining their preferred workflow, regardless of location.



HIGHLIGHTS


- Near instant streaming of images to any viewing location
- Multi-monitor application or zero-footprint browser-based viewer
- Available on Windows and iOS
- Connects to any 3rd party PACS
- Multi-vendor mammography support: Hologic, GE, Siemens, FujiFilm
- Multi-vendor imaging support: Ultrasound, X-ray, CT, MRI, mammography/tomosynthesis, ABUS
- Comprehensive image viewing and manipulation tools
- Customizable hanging protocols with mammography workflow
- Multiple concurrent users with role-based permissions
- Integrated reporting and inter-institutional image sharing



Image Viewing Functionality

Extensive viewing capabilities ensure high-quality image review with optimized display and manipulation tools across all imaging modalities.

- **Window/Level:** Auto and manual adjustment
- **Image Manipulation:** Pan/zoom, rotate, flip, invert, quadrant zoom
- **Display Features:** DICOM overlays, show/hide text, double-click min-max, MiniMap
- **Multi-vendor Support:** Including mammo/tomo studies from Hologic, GE, Siemens, FujiFilm
- **Advanced Features:** Auto-scale/align, GSPS, Secondary Capture, Mammo CAD
- **Multi-monitor support**



Professional Tools

Complete measurement and annotation toolkit enhances accuracy with efficient workflow tools and customizable interface options.

- **Measurements:** Length, angle, area (ellipse/free-hand), Cobb angles
- **Annotations:** Arrow, text, line tools with radial menu access
- **Navigation:** Magnifying glass, continuous zoom, 3D cursor, reference lines
- **Workflow:** Cine tool with adjustable speed, zoom/pan synchronization
- **Customization:** Tool layout, keypad support and customizations

Hanging Protocols

Intelligent protocol management streamlines radiologist workflow with automated layouts and customizable configurations for optimal reading efficiency.

- **Dynamic Viewport Management:** Multiple viewport setup, tiled windows, drag/drop images
- **Protocol Editor:** Customizable hanging protocols, loading/saving capabilities, single-click next-steps
- **Mammography Workflow:** Specialized designer, view modifier support
- **Intelligent Protocol Matching:** Automated layout based on study type

User Management & Collaboration

Multi-user environment with personalized settings and secure image sharing capabilities enables seamless collaboration across healthcare networks.

- **Multi-User:** Multiple concurrent users with role-based settings and permissions
- **Personalization:** Individual user preferences, roaming settings
- **Reporting:** Integrated reports with templates, key image import
- **Image Sharing:** Inter-institutional sharing, referring physician notifications, CD export

