

IHE Integration Statement			
Vendor	Product Name	Version	Date
Candelis, Inc.	Candelis ImageGrid	3.1.0	1/23/18

This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below :

Integration Profile	Actor	Integration Profile Option
Access to Radiology Information	Image Manager/Archive	None
Audit Trail and Node Authentication	Secure Node	None
Audit Trail and Node Authentication	Secure Node	Radiology Audit Trail
Consistent Presentation of Images	Image Manager/Archive	None
Consistent Time	Time Client	None
Chest X-Ray CAD Display	Image Manager/Archive	None
Digital Breast Tomosynthesis	Image Manager/Archive	None
Digital Breast Tomosynthesis	Image Manager/Archive	Key Images
Digital Breast Tomosynthesis	Image Manager/Archive	User Annotation
Digital Breast Tomosynthesis	Image Manager/Archive	For Presentation Breast Projection X-Ray Images
Digital Breast Tomosynthesis	Image Manager/Archive	For Processing Breast Projection X-Ray Images
MR Diffusion Imaging	Image Manager/Archive	None
Evidence Documents	Image Manager/Archive	None
Import Reconciliation Workflow	Image Manager/Archive	None

Internet address for vendor's IHE information :
[IHE Integration Statement 2018](#)

Links to Standards Conformance Statements for the Implementation

HL7	HL7 Conformance Statement - Redwood
DICOM	DICOM Conformance Statement - Redwood

Link to general information on IHE: <http://www.ihe.net>

IHE Integration Statement			
Vendor	Product Name	Version	Date
Candelis, Inc.	Candelis ImageGrid	3.1.0	1/23/18

This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below :

Integration Profile	Actor	Integration Profile Option
Key Image Note	Image Manager/Archive	None
Mammography Image	Image Manager/Archive	None
NM Image	Image Manager/Archive	None
Portable Data for Imaging	Portable Media Importer	DVD Media
Portable Data for Imaging	Portable Media Importer	None
Stereotactic Mammography Image	Image Manager/Archive	None
Scheduled Workflow.b	Image Manager/Archive	None
Scheduled Workflow.b	Image Manager/Archive	Performed Work Status Update Receive
Scheduled Workflow.b	Image Manager/Archive	PPS Exception Management
Scheduled Workflow.b	Performed Procedure Step Manager	None
Cross-Enterprise Document Sharing	Document Registry	None
Cross-Enterprise Document Sharing	Document Registry	Patient Identity Feed (HL7 V2)
Cross-Enterprise Document Sharing	Document Repository	None

Internet address for vendor's IHE information :
[IHE Integration Statement 2018](#)

Links to Standards Conformance Statements for the Implementation

HL7	HL7 Conformance Statement - Redwood
DICOM	DICOM Conformance Statement - Redwood

Link to general information on IHE: <http://www.ihe.net>